



Installation and operating instructions

AUL 04._L – Stainless steel trough

BASIC TECHNICAL INFORMATION

<u>Inner dimensions of the sink:</u>	AUL 04.1L	1250 (1190) x 430 (330) x 420 (188÷200)
	AUL 04.2L	1900 (1840) x 430 (330) x 420 (180÷200)
	AUL 04.3L	2500 (2440) x 430 (330) x 420 (175÷200)
	AUL 04.4L	2995 (2940) x 430 (330) x 420 (170÷200)
<u>Diameter of the outlet:</u>	D 40 mm	
<u>Weblink:</u>	AUL 04._L	


Description of the stainless steel trough

- AUL 04._L is a stainless steel wall-mounted washing trough.
- Its construction enables installation of outflow arms, or sensor-controlled taps. The standard version of the trough is without openings for the taps.
- Version type of the trough is left or right – as per the location of water outlet (from a front view).
- Medical trough is made of stainless steel material AISI 316. Trough is cased from the front and on both sides.

Installation of the stainless steel trough

Requirements for setting up construction

- Set-up the outlet (D) according to technical draw.



Assembly of the trough

1. Set up the carrying bar on the wall with anchoring bolts into the joggles. Bolts and joggles supplied by producer enable fixing into the sufficiently bearing masonry wall (solid brick or concrete wall, etc.). For fixing into other type of walls (cavity brick wall, etc.) it is necessary to use suitable anchoring material.
2. Fix the trough properly on the carrying bar. Check the horizontal level on the upper edge.



3. Fix both sides of the trough with anchoring bolts into the joggles.
4. Set-up siphon.
5. Seal the gap between the trough and the wall with suitable gap-filling glue to prevent water leaking.

Complete delivery

- Stainless steel trough 1 pc.
- Anchoring bolts, joggles and washers 7 - 10 pcs
- Plasticsiphon 1 pcs

Maintenance and cleaning

- The device was made from the stainless steel of quality corresponding with ČSN 17 346 (AISI 316) standard and, therefore, it must not be operated in chemically-aggressive environment and
- **Preparations containing chlorine must not be used for its cleaning!!**
- Cleaning agents by WÜRTH are recommended:
 - Metal renewal agent - Order No. 893 121 1
 - Stainless steel spray treatment - Order No. 0893 121 – K.

Note

- When wash trough is used in the same surrounding in which standard steel is processed, superficial corrosion can appear due to the sediments of metal particles on the bottom of the trough. It is necessary to rinse well the trough after every usage to remove this possible deposits of corrosion.